

# Calcyphosine 2 antibody



Catalog Number: 108820

## Product name

Calcyphosine 2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

Calcyphosine 2 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 54 kDa

## Preparation

This antibody was obtained by immunization of Calcyphosine 2 recombinant protein (Accession Number: XM\_005269194). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108820(CAPS2 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:108820(CAPS2 antibody) at dilution of 1:50 (under 10x lens)