

# CALML6 antibody

Catalog Number: 108829

## Product name

CALML6 antibody

## Specificity

Human; other species not tested.

## Antibody description

CALML6 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue. Positive IHC detected in human heart tissue, human prostate hyperplasia tissue. Observed molecular weight by Western-blot: 20-25 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_138705). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:108829(CALML6 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human heart using Catalog No:108829(CALML6 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:108829(CALML6 antibody) at dilution of 1:50 (under 40x lens)