

# Caspase 8/p18 antibody



Catalog Number: 108879

## General Information

### Product name:

Caspase 8/p18 antibody

### Immunogen:

[Human Caspase 8 Recombinant protein \(GST tag\)](#)

### Applications:

ELISA, WB, IHC, IP, IF

### Specificity:

Human, Mouse; other species not tested.

### Antibody description:

Caspase 8/p18 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human lymphoma tissue, human colon tissue. Positive IP detected in HeLa cells. Positive WB detected in HepG2 cells, HeLa cells, Raji cells, RAW264.7 cells, UV treated Hela cells. Observed molecular weight by Western-blot: 50-55 kDa, 25-32 kDa

## Properties

### Preparation:

This antibody was obtained by immunization of Caspase 8/p18 recombinant protein (Accession Number: BC028223). Purification method: Antigen affinity purified.

### Formulation:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

### Storage:

Store at -20°C. DO NOT ALIQUOT

### Clonality:

Polyclonal

**Ig Type:**Rabbit IgG

## Applications

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000


IHC: 1:20-1:200

IF: 1:10-1:100

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:108879(CASP8 antibody) at dilution of 1:1000




IP Result of anti-CASP8 (IP:Catalog No:108879, 3ug; Detection:Catalog No:108879 1:500) with HeLa cells lysate 2500ug.


# Caspase 8/p18 antibody

Catalog Number: 108879


---



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:108879(CASP8 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells using Catalog No:108879(Caspase 8 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:108879(CASP8 antibody) at dilution of 1:50 (under 10x lens)