

Catalase-Specific antibody



Catalog Number: 108885

Product name

Catalase-Specific antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

Human; other species not tested.

Validations

Antibody description

Catalase-Specific Rabbit Polyclonal antibody.
Positive WB detected in HT-1080 cells, HeLa cells.
Positive IHC detected in human hepat°Cirrhosis tissue. Positive IF detected in Hela cells, HepG2 cells. Observed molecular weight by Western-blot: 60kd




Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001752).
Purification method: Antigen affinity purified.

HT-1080 cells were subjected to SDS PAGE followed by western blot with Catalog No:108885(CAT-Specific antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IHC

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000


Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:108885(CAT-Specific antibody) at dilution of 1:100 (under 10x lens)



Catalase-Specific antibody

Catalog Number: 108885

Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:108885(CAT-Specific antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HeLa cells, using CAT antibody Catalog No:108885 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).