

Cathepsin A antibody



Catalog Number: 108886

Product name

Cathepsin A antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Cathepsin A Rabbit Polyclonal antibody. Positive WB detected in Sp2/0 cells. Observed molecular weight by Western-blot: 54-60 kDa

Preparation

This antibody was obtained by immunization of Cathepsin A recombinant protein (Accession Number: NM_001127695). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



Sp2/0 cells were subjected to SDS PAGE followed by western blot with Catalog No:108886(CTSA antibody) at dilution of 1:300