

Cathepsin V antibody



Catalog Number: 108894

Product name

Cathepsin V antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Cathepsin V Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, Jurkat cells, MCF-7 cells, mouse testis tissue. Observed molecular weight by Western-blot: 37-39kd, 30kd

Preparation

This antibody was obtained by immunization of Cathepsin V recombinant protein (Accession Number: NM_001333). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108894(CTSL2 antibody) at dilution of 1:800