

CCDC40 antibody



Catalog Number: 108973

Product name

CCDC40 antibody

Specificity

Human; other species not tested.

Antibody description

CCDC40 Rabbit Polyclonal antibody. Positive WB detected in fetal human brain tissue, A549 cells. Observed molecular weight by Western-blot: 90-100 kDa

Preparation

This antibody was obtained by immunization of CCDC40 recombinant protein (Accession Number: NM_001243342). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108973(CCDC40 Antibody) at dilution of 1:1000