

CCDC70 antibody



Catalog Number: 108986

Product name

CCDC70 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CCDC70 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 29 kDa

Preparation

This antibody was obtained by immunization of CCDC70 recombinant protein (Accession Number: BC069770). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:108986(CCDC70 antibody) at dilution of 1:500



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108986(CCDC70 Antibody) at dilution of 1:2500