

CCDC83 antibody



Catalog Number: 108989

Product name

CCDC83 antibody

Specificity

Human; other species not tested.

Antibody description

CCDC83 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 49 kDa

Preparation

This antibody was obtained by immunization of CCDC83 recombinant protein (Accession Number: NM_001286159). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:108989(CCDC83 Antibody) at dilution of 1:600