

CD32 antibody



Catalog Number: 109026

Product name

CD32 antibody

Specificity

Human; other species not tested.

Antibody description

CD32 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, U-937 cells. Observed molecular weight by Western-blot: 45-49 kDa

Preparation

This antibody was obtained by immunization of CD32 recombinant protein (Accession Number: NM_021642). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109026(FCGR2A antibody) at dilution of 1:500