

# CCT8 antibody



Catalog Number: 109037

## Product name

CCT8 antibody

## Specificity

Human, Mouse, Rat, Zebrafish; other species not tested.

## Antibody description

CCT8 Rabbit Polyclonal antibody. Positive IF detected in MCF-7 cells. Positive WB detected in mouse kidney tissue, Jurkat cells, Raji cells. Observed molecular weight by Western-blot: 60kd

## Preparation

This antibody was obtained by immunization of CCT8 recombinant protein (Accession Number: NM\_001282908). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IF, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



Immunofluorescent analysis of MCF-7 cells, using CCT8 antibody Catalog No:109037 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:109037(CCT8 Antibody) at dilution of 1:1500