

CD11B/Integrin alpha M antibody



Catalog Number: 109042

Product name

CD11B/Integrin alpha M antibody

Specificity

Human; other species not tested.

Antibody description

CD11B/Integrin alpha M Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue, human spleen tissue.

Preparation

This antibody was obtained by immunization of CD11B/Integrin alpha M recombinant protein (Accession Number: NM_000632). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109042(CD11B Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109042(CD11B Antibody) at dilution of 1:50 (under 40x lens)