

CCR5 antibody

Catalog Number: 109077

Product name

CCR5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CCR5 Rabbit Polyclonal antibody. Positive FC detected in K-562 cells. Positive WB detected in HeLa cells, K-562 cells. Observed molecular weight by Western-blot: 45 kDa

Preparation

This antibody was obtained by immunization of CCR5 recombinant protein (Accession Number: NM_001100168). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, FC

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109077(CCR5 antibody) at dilution of 1:800

Plot P03, gated on P01.R1



1X10⁶ K-562 cells were stained with 0.2ug CCR5 antibody (Catalog No:109077, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.