

CDC14A antibody



Catalog Number: 109091

Product name

CDC14A antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

CDC14A Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, HeLa cells, human brain tissue. Observed molecular weight by Western-blot: 66kd

Preparation

This antibody was obtained by immunization of CDC14A recombinant protein (Accession Number: BC038979). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:109091(CDC14A antibody) at dilution of 1:500