

Cdc20 antibody

Catalog Number: 109094

Product name

Cdc20 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Cdc20 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HEK-293 cells, HL-60 cells, PC-3 cells. Positive IP detected in Jurkat cells. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of Cdc20 recombinant protein (Accession Number: NM_001255). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109094(CDC20 Antibody) at dilution of 1:800



IP Result of anti-CDC20 (IP:Catalog No:109094, 4ug; Detection:Catalog No:109094 1:800) with Jurkat cells lysate 3200ug.