

CDC42SE2 antibody



Catalog Number: 109108

Product name

CDC42SE2 antibody

Specificity

Human; other species not tested.

Antibody description

CDC42SE2 Rabbit Polyclonal antibody. Positive IHC detected in human colon tissue, human lymphoma tissue.

Preparation

This antibody was obtained by immunization of CDC42SE2 recombinant protein (Accession Number: XM_047417398). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:109108(CDC42SE2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:109108(CDC42SE2 Antibody) at dilution of 1:50 (under 40x lens)