

# CDCA8 antibody

Catalog Number: 109118

## Product name

CDCA8 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CDCA8 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue. Positive WB detected in mouse kidney tissue. Observed molecular weight by Western-blot: 38-40 kDa

## Preparation

This antibody was obtained by immunization of CDCA8 recombinant protein (Accession Number: NM\_001256875). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:109118(CDCA8 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:109118(CDCA8 antibody) at dilution of 1:50 (under 10x lens)