

CD58 antibody



Catalog Number: 109130

Product name

CD58 antibody

Specificity

Human; other species not tested.

Antibody description

CD58 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human lymph tissue, Jurkat cells. Observed molecular weight by Western-blot: 65-70 kDa

Preparation

This antibody was obtained by immunization of CD58 recombinant protein (Accession Number: BC005930). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109130(CD58 antibody) at dilution of 1:500