

CDK1 antibody

Catalog Number: 109155

Product name

CDK1 antibody

Specificity

Human; other species not tested.

Antibody description

CDK1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 34kd, 27kd

Preparation

This antibody was obtained by immunization of CDK1 recombinant protein (Accession Number: NM_001786). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109155(CDK1 antibody) at dilution of 1:300