

CDK3 antibody



Catalog Number: 109160

Product name

CDK3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CDK3 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, Jurkat cells, mouse brain tissue. Observed molecular weight by Western-blot: 40-45kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001258). Purification method: Antigen affinity purified.

Formulation

0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:109160(CDK3 antibody) at dilution of 1:800