

# CDK6 antibody

Catalog Number: 109166

## Product name

CDK6 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CDK6 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells, C6 cells, K-562 cells. Observed molecular weight by Western-blot: 37-40 kDa

## Preparation

This antibody was obtained by immunization of CDK6 recombinant protein (Accession Number: BC027989). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109166(CDK6 antibody) at dilution of 1:2000



IP Result of anti-CDK6 (IP:Catalog No:109166, 4ug; Detection:Catalog No:109166 1:700) with Jurkat cells lysate 3200ug.