

CEP78 antibody



Catalog Number: 109187

Product name

CEP78 antibody

Specificity

Human; other species not tested.

Antibody description

CEP78 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_032171).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:50-1:500

Validations



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using Catalog No:109187(CEP78 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)