

Ceruloplasmin antibody



Catalog Number: 109191

Product name

Ceruloplasmin antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Ceruloplasmin Rabbit Polyclonal antibody. Positive WB detected in human blood tissue, mouse brain tissue. Observed molecular weight by Western-blot: 140kd,85kd

Preparation

This antibody was obtained by immunization of Ceruloplasmin recombinant protein (Accession Number: NM_000096). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human blood tissue were subjected to SDS PAGE followed by western blot with Catalog No:109191(CP antibody) at dilution of 1:1000