

CETP antibody

Catalog Number: 109196

Product name

CETP antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CETP Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human brain tissue, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 66-75kd

Preparation

This antibody was obtained by immunization of CETP recombinant protein (Accession Number: NM_000078). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications



ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



WB result of Catalog No:109196.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109196(CETP antibody) at dilution of 1:600