

# CHERP antibody



Catalog Number: 109223

## Product name

CHERP antibody

## Specificity

Human; other species not tested.

## Antibody description

CHERP Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Positive IHC detected in mouse brain tissue. Observed molecular weight by Western-blot: 130 kDa

## Preparation

This antibody was obtained by immunization of CHERP recombinant protein (Accession Number: BC021294). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109223(CHERP Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:109223(CHERP Antibody) at dilution of 1:50 (under 40x lens)