

# CHPF-N terminal antibody



Catalog Number: 109274

## Product name

CHPF-N terminal antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CHPF-N terminal Rabbit Polyclonal antibody.  
Positive WB detected in human brain tissue.  
Observed molecular weight by Western-blot: 85-90 kDa, 65-70 kDa

## Preparation

This antibody was obtained by immunization of CHPF-N terminal recombinant protein (Accession Number: BC008878). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109274(CHPF antibody) at dilution of 1:800