

Chromogranin A antibody



Catalog Number: 109287

Product name

Chromogranin A antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Chromogranin A Rabbit Polyclonal antibody. Positive IHC detected in human pancreas tissue, human colon tissue, human lung cancer tissue, mouse colon tissue.

Preparation

This antibody was obtained by immunization of Chromogranin A recombinant protein (Accession Number: NM_001275). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500


Validations



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using Catalog No:109287(CHGA Antibody) at dilution of 1:400 (under 10x lens).




Immunohistochemistry of paraffin-embedded human pancreas tissue slide using Catalog No:109287(CHGA Antibody) at dilution of 1:400 (under 40x lens).



Immunohistochemical of paraffin-embedded human colon using Catalog No:109287(CHGA antibody) at dilution of 1:50 (under 10x lens)

Chromogranin A antibody

Catalog Number: 109287



Immunohistochemical of paraffin-embedded human colon using Catalog No:109287(CHGA antibody) at dilution of 1:50 (under 40x lens)