

CHST10 antibody



Catalog Number: 109291

Product name

CHST10 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CHST10 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive WB detected in BxPC-3 cells. Observed molecular weight by Western-blot: 57 kDa

Preparation

This antibody was obtained by immunization of CHST10 recombinant protein (Accession Number: XM_017005381). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109291(CHST10 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:109291(CHST10 antibody) at dilution of 1:50 (under 10x lens)