

CHST13-Specific antibody



Catalog Number: 109294

Product name

CHST13-Specific antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Antibody description

CHST13-Specific Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive IHC detected in human kidney tissue, human hepat°Cirrhosis tissue. Positive IP detected in mouse liver tissue. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 25 kDa

Validations



Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_152889). Purification method: Antigen affinity purified.

mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109294(CHST13-Specific antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IF, IP

Dilutions

Recommended Dilution:


WB: 1:500-1:5000




Immunohistochemical of paraffin-embedded human kidney using Catalog No:109294(CHST13-Specific antibody) at dilution of 1:100 (under 10x lens)

CHST13-Specific antibody

Catalog Number: 109294




Immunohistochemical of paraffin-embedded human kidney using Catalog No:109294(CHST13-Specific antibody) at dilution of 1:100 (under 40x



lens)

Immunofluorescent analysis of HeLa cells, using CHST13 antibody Catalog No:109294 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-CHST13-Specific (IP:Catalog No:109294, 4ug; Detection:Catalog No:109294 1:500) with mouse liver tissue lysate 4000ug.