

CHSY1 antibody



Catalog Number: 109301

Product name

CHSY1 antibody

Specificity

human, mouse, rat, zebrafish; other species not tested.

Antibody description

CHSY1 Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue, mouse heart tissue. Positive IHC detected in human lung cancer. Observed molecular weight by Western-blot: 66kd

Preparation

This antibody was obtained by immunization of CHSY1 recombinant protein (Accession Number: NM_014918). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:109301(CHSY1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:109301(CHSY1 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:109301(CHSY1 antibody) at dilution of 1:100

CHSY1 antibody



Catalog Number: 109301

(under 40x lens)