

# CIN85 antibody



Catalog Number: 109316

## Product name

CIN85 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CIN85 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, fetal human brain tissue. Observed molecular weight by Western-blot: 85 kDa

## Preparation

This antibody was obtained by immunization of CIN85 recombinant protein (Accession Number: NM\_031892). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109316(CIN85 Antibody) at dilution of 1:1000



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109316(CIN85 Antibody) at dilution of 1:1000