

CITED2 antibody



Catalog Number: 109324

Product name

CITED2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CITED2 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, mouse brain tissue. Positive IF detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 24 kDa

Preparation

This antibody was obtained by immunization of CITED2 recombinant protein (Accession Number: NM_001168388). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:50-1:500

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:109324(CITED2 antibody) at dilution of 1:400



Immunofluorescent analysis of (10% Formaldehyde) fixed SH-SY5Y cells using Catalog No:109324(CITED2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)