

# Claudin 23 antibody



Catalog Number: 109340

## Product name

Claudin 23 antibody

## Specificity

Human; other species not tested.

## Antibody description

Claudin 23 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 32 kDa

## Preparation

This antibody was obtained by immunization of Claudin 23 recombinant protein (Accession Number: NM\_194284). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:109340(CLDN23 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:109340(CLDN23 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:109340(CLDN23 Antibody) at dilution of 1:50 (under 40x lens)