

CNTN1 antibody

Catalog Number: 109362

Product name

CNTN1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CNTN1 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue. Positive IP detected in rat brain tissue. Observed molecular weight by Western-blot: 135 kDa

Preparation

This antibody was obtained by immunization of CNTN1 recombinant protein (Accession Number: XM_017018825). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109362(CNTN1 Antibody) at dilution of 1:1200



IP Result of anti-CNTN1 (IP:Catalog No:109362, 4ug; Detection:Catalog No:109362 1:1500) with rat brain tissue lysate 4000ug.