

Connexin-26 antibody



Catalog Number: 109449

Product name

Connexin-26 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Connexin-26 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, mouse brain tissue. Observed molecular weight by Western-blot: 26kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004004). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:109449(Connexin-26 antibody) at dilution of 1:500