

COPZ1 antibody

Catalog Number: 109467

Product name

COPZ1 antibody

Specificity

Human, Mouse, Rat, Monkey; other species not tested.

Antibody description

COPZ1 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, COS-7 cells, Jurkat cells. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 20 kDa

Preparation

This antibody was obtained by immunization of COPZ1 recombinant protein (Accession Number: NM_016057). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:50-1:500

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:109467(COPZ1 antibody) at dilution of 1:500



Immunofluorescent analysis of (-20°C Ethanol) fixed HEK-293 cells using Catalog No:109467(COPZ1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)