

COQ5 antibody



Catalog Number: 109470

Product name

COQ5 antibody

Immunogen

[Human COQ5 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

COQ5 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, mouse brain tissue, mouse liver tissue. Positive IHC detected in human heart tissue, human brain tissue, human kidney tissue, human liver tissue, human lung tissue, human ovary tissue, human placenta tissue, human skin tissue, human spleen tissue, human testis tissue. Observed molecular weight by Western-blot: 34.7kd

Preparation

This antibody was obtained by immunization of COQ5 recombinant protein (Accession Number: NM_032314). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations




mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:109470(COQ5 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human heart using Catalog No:109470(COQ5 antibody) at dilution of 1:100 (under 10x lens)

COQ5 antibody

Catalog Number: 109470



Immunohistochemical of paraffin-embedded human heart using Catalog No:109470(COQ5 antibody) at dilution of 1:100 (under 40x lens)