

# CPNE5 antibody



Catalog Number: 109517

## Product name

CPNE5 antibody

## Immunogen

[Human CPNE5 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CPNE5 Rabbit Polyclonal antibody. Positive WB detected in mouse embryo tissue. Observed molecular weight by Western-blot: 68 kDa

## Preparation

This antibody was obtained by immunization of CPNE5 recombinant protein (Accession Number: NM\_001376890). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



mouse embryo tissue were subjected to SDS PAGE followed by western blot with Catalog No:109517(CPNE5 Antibody) at dilution of 1:300