

Creatine Kinase antibody



Catalog Number: 109541

Product name

Creatine Kinase antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Creatine Kinase Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 43kd

Preparation

This antibody was obtained by immunization of Creatine Kinase recombinant protein (Accession Number: NM_001824). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:109541 (CKM antibody) at dilution of 1:1000