

CST9L antibody



Catalog Number: 109615

Product name

CST9L antibody

Specificity

human; other species not tested.

Antibody description

CST9L Rabbit Polyclonal antibody. Positive WB detected in human kidney tissue. Observed molecular weight by Western-blot: 17kd

Preparation

This antibody was obtained by immunization of CST9L recombinant protein (Accession Number: NM_080610). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:109615(CST9L antibody) at dilution of 1:500