

CSTL1-Specific antibody



Catalog Number: 109620

Product name

CSTL1-Specific antibody

Specificity

human, mouse; other species not tested.

Antibody description

CSTL1-Specific Rabbit Polyclonal antibody.
Positive IHC detected in human testis.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_138283).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human testis using Catalog No:109620(CSTL1-Specific antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human testis using Catalog No:109620(CSTL1-Specific antibody) at dilution of 1:50 (under 40x lens)