

CWF19L2 antibody



Catalog Number: 109632

Product name

CWF19L2 antibody

Immunogen

[Human CWF19L2 Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

CWF19L2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, Jurkat cells, PC-3 cells. Positive IHC detected in human kidney tissue, human brain tissue, human breast cancer tissue. Observed molecular weight by Western-blot: 103 kDa

Preparation

This antibody was obtained by immunization of CWF19L2 recombinant protein (Accession Number: BC110442). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations




HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109632(CWF19L2 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:109632(CWF19L2 Antibody) at dilution of 1:50 (under 40x lens)

CWF19L2 antibody

Catalog Number: 109632



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:109632(CWF19L2 Antibody) at dilution of 1:50 (under 40x lens)