

# CX3CL1 antibody

Catalog Number: 109633

## General Information

### Product name:

CX3CL1 antibody

### Immunogen:

[Human CX3CL1 Recombinant protein \(GST tag\)](#)

### Applications:

ELISA, WB, IHC

### Specificity:

Human; other species not tested.

### Antibody description:

CX3CL1 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, human brain tissue. Positive IHC detected in human small intestine tissue. Observed molecular weight by Western-blot: 90-100 kDa

## Properties

### Preparation:

This antibody was obtained by immunization of CX3CL1 recombinant protein (Accession Number: BC001163). Purification method: Antigen affinity purified.

### Formulation:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

### Storage:

Store at -20°C. DO NOT ALIQUOT

### Clonality:

Polyclonal

**Ig Type:**Rabbit IgG

## Applications

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109633(CX3CL1 antibody) at dilution of 1:1000




Immunohistochemistry of paraffin-embedded human small intestine slide using Catalog No:109633(CX3CL1 Antibody) at dilution of 1:50 (under 10x lens)

# CX3CL1 antibody

Catalog Number: 109633

---



Immunohistochemistry of paraffin-embedded human small intestine slide using Catalog No:109633(CX3CL1 Antibody) at dilution of 1:50 (under 40x lens)