

CXCR4 antibody

Catalog Number: 109644

Product name

CXCR4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CXCR4 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in rat spleen tissue, HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 60-70kd

Preparation

This antibody was obtained by immunization of CXCR4 recombinant protein (Accession Number: NM_003467). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



IP Result of anti-CXCR4 (IP:Catalog No:109644, 3ug; Detection:Catalog No:109644 1:500) with



HeLa cells lysate 1000ug.

rat spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:109644(CXCR4 Antibody) at dilution of 1:1000