

CYC1 antibody

Catalog Number: 109659

Product name

CYC1 antibody

Specificity

Human; other species not tested.

Antibody description

CYC1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 35 kDa

Preparation

This antibody was obtained by immunization of CYC1 recombinant protein (Accession Number: NM_001916). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109659(CYC1 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:109659(CYC1 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:109659(CYC1 antibody) at dilution of 1:100 (under 40x lens)