

CYP4A11 antibody



Catalog Number: 109703

Product name

CYP4A11 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

CYP4A11 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, rat kidney tissue. Positive IHC detected in human kidney. Observed molecular weight by Western-blot: 59kd

Preparation

This antibody was obtained by immunization of CYP4A11 recombinant protein (Accession Number: BC022851). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:109703(CYP4A11 antibody) at dilution of 1:200



Immunohistochemical of paraffin-embedded human kidney using Catalog No:109703(CYP4A11 antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:109703(CYP4A11 antibody) at dilution of 1:50 (under 40x lens)