

# DCLK1 antibody



Catalog Number: 109751

## Product name

DCLK1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DCLK1 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue, human brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 48 kDa, 82 kDa

## Preparation

This antibody was obtained by immunization of DCLK1 recombinant protein (Accession Number: NM\_004734). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109751(DCLK1 Antibody) at dilution of 1:1500