

DCLRE1C antibody



Catalog Number: 109753

Product name

DCLRE1C antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DCLRE1C Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human kidney tissue. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of DCLRE1C recombinant protein (Accession Number: BC022254). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109753(DCLRE1C antibody) at dilution of 1:400