

# DCLRE1C antibody



Catalog Number: 109753

## Product name

DCLRE1C antibody

Polyclonal

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat; other species not tested.

## Applications

ELISA, WB

## Antibody description

DCLRE1C Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human kidney tissue. Observed molecular weight by Western-blot: 60 kDa

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of DCLRE1C recombinant protein (Accession Number: BC022254). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109753(DCLRE1C antibody) at dilution of 1:400

## Clonality