

# DCUN1D1 antibody



Catalog Number: 109763

## Product name

DCUN1D1 antibody

ELISA, IHC, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

IHC: 1:20-1:200


## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DCUN1D1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, mouse kidney tissue. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 28 kDa

## Validations




## Preparation

This antibody was obtained by immunization of DCUN1D1 recombinant protein (Accession Number: NM\_020640). Purification method: Antigen affinity purified.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109763(DCUN1D1 antibody) at dilution of

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



1:300

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

Immunohistochemical of paraffin-embedded human gliomas using Catalog No:109763(DCUN1D1 antibody) at dilution of 1:50 (under 10x lens)

## Applications